

Recently, one of our clients upgraded to Meditech 6.07 PP12 and we discovered that Meditech has now added Event Types in the Outbound ADT Interface. For anyone who has worked with Meditech interfaces, you know this is a welcomed feature. As we worked with Meditech, we were able to ascertain a list of Event Types that they have included in the interface. These values will populate EVN.4.



There are many ways that this can be utilize in your ADT interfaces. The first that comes to mind is filtering specific events. For example, we have run across many times when working with Care Fusion/MedMined that they only want A08 message that include Height/Weight changes. In EVN.4, you will find NMIHGTWGT for Enter/Edit Height and Weight. You can filter on MSH.9 = ADT^A08 and EVN.4 = NMIHGTWGHT and now send just the message that contain a height and weight enter edit.

The following is a list of Meditech values that can appear in EVN.4:

- "EDACCT#" - "EDIT ACCOUNT #"
 - "EDUNIT#" - "EDIT UNIT #"
 - "EDSTATUS" - "EDIT PRE/SCH STATUS"
 - "MERGE" - "MRI MERGE"
 - "UNMERGE" - "MRI UNMERGE"
 - "EEVIPCON" - "EDIT VIP/CONFIDENTIAL STATUS" ,
 - "ADDOUTLOC" - "ADD "_((.S).PTO@Tr.l2u)" LOCATION" ,
 - "DELOUTLOC" - "DELETE "_((.S).PTO@Tr.l2u)" LOCATION" ,
 - "EDPRIMLOC" - "EDIT PRIMARY LOCATION" ,
 - "EDGUAR" - "EDIT GUARANTOR DATA" ,
 - "EDINS" - "EDIT INSURANCE DATA" ,
 - "EDNAM" - "EDIT NAME" ,
 - "EDSEX" - "EDIT SEX" ,
 - "EDAGE" - "EDIT AGE" ,
 - "EDBD" - "EDIT BIRTHDATE" ,
 - "EDLOC" - "EDIT LOCATION" ,
 - "EDZIP" - "EDIT ZIP/POSTAL CODE" ,
 - "FXPCDR" - "EDIT PRIMARY CARE DOCTOR"
 - "PREDISREF" - "PRE DISCHARGE REFERRAL"
-
- "DEL" - "DELETE "
 - "ABS" - "ABSTRACTING DATA" ,
 - "NSC" - "NURSING DATA"
 - "NSC2" - "NURSING DATA (CONT)"
 - "NSCTL" - "NURSING DATA TEMP LOC"
 - "NMIAL1" - "EDIT ALLERGY"

"NMIEN2" - "ENTER ADM/REG DATA (CONT)"
"NMIED2" - "EDIT ADM/REG DATA (CONT)"
"NMIGUAR" - "GUARANTOR DATA"
"NMIGEMP" - "GUAR EMPLOYER DATA"
"NMIHGTWGT" - "ENTER/EDIT HEIGHT & WEIGHT"
"NMIINS" - "INSURANCE DATA"
"NMIINS2" - "INSURANCE DATA (CONT)"
"NMIINSQ" - "INSURANCE QUERY DATA"
"NMIINSU" - "INSURANCE UB92 DATA"
"NMIUB82" - "UB82 DATA"
"NMIUB822" - "UB82 DATA (CONT)"

Inpatients - SCH and PRE have no event.stats.program

"ENSCHIN" - "ENTER SCH IN"
"EDSCHIN" - "EDIT SCH IN"
"ENSCHOB" - "ENTER SCH OBSERV"
"EDSCHOB" - "EDIT SCH OBSERV"
"ENPREIN" - "ENTER PRE IN"
"EDPREIN" - "EDIT PRE IN"
"ENPREOB" - "ENTER PRE OBSERV"
"EDPREOB" - "EDIT PRE OBSERV"
"ENPRENB" - "ENTER PRE NEWBORN"
"EDPRENB" - "EDIT PRE NEWBORN"
"ENADMIN" - "ENTER ADM IN"
"ENADMOB" - "ENTER OBSERV"
"ENADMNB" - "ENTER NEWBORN"
"EDADMIN" - "EDIT ADM IN"
"ENDISIN" - "DISCHARGE IN"
"ENDISOB" - "DISCHARGE OBSERV"
"EDDISIN" - 004^@event.oth.vendor.tx.type,"IN"^@event.pt.type,
"EDIT DIS IN" ,"eddisin"^@event.stats.program,
"FXADMDR" - 004^@event.oth.vendor.tx.type,"IN"^@event.pt.type,
"EDIT ADMIT DOCTOR" ,"fxadmdr"^@event.stats.program,
"FXATTDR" - 004^@event.oth.vendor.tx.type,"IN"^@event.pt.type,
"EDIT ATTEND DOCTOR" ,"fxattdr"^@event.stats.program,
"FXFC" - 004^@event.oth.vendor.tx.type,"IN"^@event.pt.type,
"EDIT FINANCIAL CLASS" ,"fxfc"^@event.stats.program,
"FXSVC" - 004^@event.oth.vendor.tx.type,"IN"^@event.pt.type,
"EDIT SERVICE" ,"fxsvc"^@event.stats.program,
"FXACCOM" - 004^@event.oth.vendor.tx.type,"IN"^@event.pt.type,
"EDIT REQUESTED ACCOM" ,"fxaccom"^@event.stats.program,
"FXLOA" - 004^@event.oth.vendor.tx.type,"IN"^@event.pt.type,

"EDIT LOA" , "fxloa"^^@event.stats.program, "Y"^^@event.undoable,
"FXRRACCOM" - 004^^@event.oth.vendor.tx.type, "IN"^^@event.pt.type,
"EDIT ROOM RATE ACCOM" , "fxrracom"^^@event.stats.program,
"FXDISDT" - 003^^@event.oth.vendor.tx.type, "IN"^^@event.pt.type,
"EDIT DISCHARGE DATE" , "fxdisdt"^^@event.stats.program,
"SWAPADMIN" - 007^^@event.oth.vendor.tx.type, "IN"^^@event.pt.type,
"SWAP BEDS" , "tfradmin"^^@event.stats.program, "Y"^^@event.undoable,
"TFRADMIN" - 002^^@event.oth.vendor.tx.type, "IN"^^@event.pt.type,
"CHANGE ROOM/BED" , "tfradmin"^^@event.stats.program, "Y"^^@event.undoable,
"UNDOSWAP" - 007^^@event.oth.vendor.tx.type, "IN"^^@event.pt.type,
"UNDO SWAP BEDS" , "undotfr"^^@event.stats.program,
"UNDOTFR" - 002^^@event.oth.vendor.tx.type, "IN"^^@event.pt.type,
"UNDO CHANGE ROOM/BED" , "undotfr"^^@event.stats.program,
"UNDOLOA" - 004^^@event.oth.vendor.tx.type,
"UNDO LOA" , "fxloa"^^@event.stats.program,
"INTOOBSV" - 004^^@event.oth.vendor.tx.type, "IN"^^@event.pt.type,
"CHANGE IN TO OBSERV" , "intoobsv"^^@event.stats.program,
"PREINTOOBSV" - "IN"^^@event.pt.type, "CHANGE PRE IN TO OBSERV" ,
"SCHINTOOBSV" - "IN"^^@event.pt.type, "CHANGE SCH IN TO OBSERV" ,
"OBSVTOIN" - 004^^@event.oth.vendor.tx.type, "IN"^^@event.pt.type,
"CHANGE OBSERV TO IN" , "obsvtoin"^^@event.stats.program, "Y"^^@event.undoable,
"PREOBSVTOIN" - "IN"^^@event.pt.type, "CHANGE PRE OBSERV TO IN" ,
"SCHOBSVTOIN" - "IN"^^@event.pt.type, "CHANGE SCH OBSERV TO IN" ,
"UNDOADMIN" - 099^^@event.oth.vendor.tx.type, "IN"^^@event.pt.type,
"UNDO ADMIT IN" , "undoadmin"^^@event.stats.program,
"UNDOADMOb" - 099^^@event.oth.vendor.tx.type, "IN"^^@event.pt.type,
"UNDO ADMIT OBSERV" , "undoadmob"^^@event.stats.program,
"UNDODISIN" - 012^^@event.oth.vendor.tx.type, "IN"^^@event.pt.type,
"UNDO DISCHARGE IN" , "undodisin"^^@event.stats.program,
"UNDODISOB" - 012^^@event.oth.vendor.tx.type, "IN"^^@event.pt.type,
"UNDO DISCHARGE OBSERV" , "undodisob"^^@event.stats.program,
"CANPREIN" - "CANCEL PRE IN" , 099^^@event.oth.vendor.tx.type, "IN"^^@event.pt.type,
"CANSCHIN" - "CANCEL SCH IN" , 099^^@event.oth.vendor.tx.type, "IN"^^@event.pt.type,

CLI outpatients

"ENSCHCLI" - "ENTER SCH CLI" , 001^^@event.oth.vendor.tx.type, "CLI"^^@event.pt.type,
"EDSCHCLI" - "EDIT SCH CLI" , 004^^@event.oth.vendor.tx.type, "CLI"^^@event.pt.type,
"ENPRECLI" - "ENTER PRE CLI" , 001^^@event.oth.vendor.tx.type, "CLI"^^@event.pt.type,
"EDPRECLI" - "EDIT PRE CLI" , 004^^@event.oth.vendor.tx.type, "CLI"^^@event.pt.type,
"ENREGCLI" - 001^^@event.oth.vendor.tx.type, "Y"^^@event.undoable, "CLI"^^@event.pt.type,
"ENTER REG CLI" , "eereg"^^@event.stats.program,
"EDREGCLI" - 004^^@event.oth.vendor.tx.type, "CLI"^^@event.pt.type,

"EDIT REG CLI" , "eereg"^@event.stats.program,
"FXREGCLI" - 006^@event.oth.vendor.tx.type,"CLI"^@event.pt.type,
"UNDO CLI REGISTRATION" , "fxreg"^@event.stats.program,
"CANPRECLI" - "CANCEL PRE CLI" ,006^@event.oth.vendor.tx.type,"CLI"^@event.pt.type,
"CANSCHCLI" - "CANCEL SCH CLI" ,006^@event.oth.vendor.tx.type,"CLI"^@event.pt.type,

Outpatient room/bed events

"ENCLIBED" - 002^@event.oth.vendor.tx.type,"CLI"^@event.pt.type,"Y"^@event.undoable,
"PLACE CLI IN ROOM/BED" , "enopbed"^@event.stats.program,
"EDCLIBED" - 002^@event.oth.vendor.tx.type,"CLI"^@event.pt.type,"Y"^@event.undoable,
"EDIT CLI ROOM/BED" , "edopbed"^@event.stats.program,
"RMCLIBED" - 002^@event.oth.vendor.tx.type,"CLI"^@event.pt.type,"Y"^@event.undoable,
"REMOVE CLI FROM ROOM/BED" , "rmopbed"^@event.stats.program,
"ENDEPCLI" - 003^@event.oth.vendor.tx.type,"CLI"^@event.pt.type,"Y"^@event.undoable,
"ENTER CLI DEPARTURE" , "endepout"^@event.stats.program,
"EDDEPCLI" - 004^@event.oth.vendor.tx.type,"CLI"^@event.pt.type,
"EDIT DEPARTED CLI" , "eereg"^@event.stats.program,
"FXDEPCLIDT" - ""^@event.oth.vendor.tx.type,"CLI"^@event.pt.type,
"EDIT CLI DEPART DATE" , "fxdepdt"^@event.stats.program,
"FXCLIRRACCOM" - ""^@event.oth.vendor.tx.type,"CLI"^@event.pt.type,"Y"^@event.undoable,
"EDIT CLI ROOM RATE ACCOM" , "fxoprracom"^@event.stats.program,
"FXCLISVC" - ""^@event.oth.vendor.tx.type,"CLI"^@event.pt.type,"Y"^@event.undoable,
"EDIT CLI SERVICE" , "fxopsvc"^@event.stats.program,
"FXCLIACCOM" - ""^@event.oth.vendor.tx.type,"CLI"^@event.pt.type,"Y"^@event.undoable,
"EDIT CLI REQUESTED ACCOM" , "fxopacom"^@event.stats.program,
"UNDOENCLIBED" - 002^@event.oth.vendor.tx.type,"CLI"^@event.pt.type,
"UNDO PLACE CLI IN ROOM/BED" , "undoenopbed"^@event.stats.program,
"UNDOEDCLIBED" - 002^@event.oth.vendor.tx.type,"CLI"^@event.pt.type,
"UNDO EDIT CLI ROOM/BED" , "undopedbed"^@event.stats.program,
"UNDOFXCLISVC" - ""^@event.oth.vendor.tx.type,"CLI"^@event.pt.type,
"UNDO EDIT CLI SERVICE" , "undofxopsvc"^@event.stats.program,
"UNDOFXCLIACCOM" - ""^@event.oth.vendor.tx.type,"CLI"^@event.pt.type,
"UNDO EDIT CLI REQUESTED ACCOM" , "undofxopacom"^@event.stats.program,
"UNDOFXCLIRRACCOM" - ""^@event.oth.vendor.tx.type,"CLI"^@event.pt.type,
"UNDO EDIT CLI ROOM RATE ACCOM" , "undofxoprracom"^@event.stats.program,
"UNDORMCLIBED" - ""^@event.oth.vendor.tx.type,"CLI"^@event.pt.type,
"UNDO REMOVE CLI FROM ROOM/BED" , "undormopbed"^@event.stats.program,
"UNDODEPCLI" - 012^@event.oth.vendor.tx.type,"CLI"^@event.pt.type,
"UNDO CLI DEPARTURE" , "undodepout"^@event.stats.program,

ER outpatients:

No SCH for ER patients

"ENPREER" - "ENTER PRE ER" ,001^@event.oth.vendor.tx.type,"ER"^@event.pt.type,
"EDPREER" - "EDIT PRE ER" ,004^@event.oth.vendor.tx.type,"ER"^@event.pt.type,
"ENREGER" - 001^@event.oth.vendor.tx.type,"Y"^@event.undoable,"ER"^@event.pt.type,
"ENTER REG ER" ,"eereg"^@event.stats.program,
"EDREGER" - 004^@event.oth.vendor.tx.type,"ER"^@event.pt.type,
"EDIT REG ER" ,"eereg"^@event.stats.program,
"ENDEPER" - "DEPART ER" ,003^@event.oth.vendor.tx.type,
"Y"^@event.undoable,"ER"^@event.pt.type,"endepout"^@event.stats.program,
"EDDEPER" - 004^@event.oth.vendor.tx.type,"ER"^@event.pt.type,
"EDIT DEPARTED ER" ,"eereg"^@event.stats.program,
"UNDODEPER" - "UNDO DEPART ER" ,012^@event.oth.vendor.tx.type,
"ER"^@event.pt.type,"undodepout"^@event.stats.program,
"FXREGER" - 006^@event.oth.vendor.tx.type,"ER"^@event.pt.type,
"UNDO ER REGISTRATION" ,"fxreg"^@event.stats.program,
"CANPREER" - "CANCEL PRE ER" ,006^@event.oth.vendor.tx.type,"ER"^@event.pt.type,

ER - Outpatient room/bed events

"ENERBED" - 002^@event.oth.vendor.tx.type,"ER"^@event.pt.type,"Y"^@event.undoable,
"PLACE ER IN ROOM/BED" ,"enopbed"^@event.stats.program,
"EDERBED" - 002^@event.oth.vendor.tx.type,"ER"^@event.pt.type,"Y"^@event.undoable,
"EDIT ER ROOM/BED" ,"edopbed"^@event.stats.program,
"RMERBED" - 002^@event.oth.vendor.tx.type,"ER"^@event.pt.type,"Y"^@event.undoable,
"REMOVE ER FROM ROOM/BED" ,"rmopbed"^@event.stats.program,
"FXERRRACCOM" - ""^@event.oth.vendor.tx.type,"ER"^@event.pt.type,"Y"^@event.undoable,
"EDIT ER ROOM RATE ACCOM" ,"fxoprracom"^@event.stats.program,
"FXERSVC" - ""^@event.oth.vendor.tx.type,"ER"^@event.pt.type,"Y"^@event.undoable,
"EDIT ER SERVICE" ,"fxopsvc"^@event.stats.program,
"FXERACCOM" - ""^@event.oth.vendor.tx.type,"ER"^@event.pt.type,"Y"^@event.undoable,
"EDIT ER REQUESTED ACCOM" ,"fxopaccom"^@event.stats.program,
"UNDOENERBED" - 002^@event.oth.vendor.tx.type,"ER"^@event.pt.type,
"UNDO PLACE ER IN ROOM/BED" ,"undoenopbed"^@event.stats.program,
"UNDOEDERBED" - 002^@event.oth.vendor.tx.type,"ER"^@event.pt.type,
"UNDO EDIT ER ROOM/BED" ,"undopedbed"^@event.stats.program,
"UNDOFXERSVC" - ""^@event.oth.vendor.tx.type,"ER"^@event.pt.type,
"UNDO EDIT ER SERVICE" ,"undofxopsvc"^@event.stats.program,
"UNDOFXERACCOM" - ""^@event.oth.vendor.tx.type,"ER"^@event.pt.type,
"UNDO EDIT ER REQUESTED ACCOM" ,"undofxopaccom"^@event.stats.program,
"UNDOFXERRRACCOM" - ""^@event.oth.vendor.tx.type,"ER"^@event.pt.type,
"UNDO EDIT ER ROOM RATE ACCOM" ,"undofxoprracom"^@event.stats.program,
"UNDORMERBED" - ""^@event.oth.vendor.tx.type,"ER"^@event.pt.type,
"UNDO REMOVE ER FROM ROOM/BED" ,"undormopbed"^@event.stats.program,

RCR outpatients

"ENSCHRRCR" - "ENTER SCH RCR" ,001^@event.oth.vendor.tx.type,"RCR"^@event.pt.type,
"EDSCHRRCR" - "EDIT SCH RCR" ,004^@event.oth.vendor.tx.type,"RCR"^@event.pt.type,
"ENPRERCR" - "ENTER PRE RCR" ,001^@event.oth.vendor.tx.type,"RCR"^@event.pt.type,
"EDPRERCR" - "EDIT PRE RCR" ,004^@event.oth.vendor.tx.type,"RCR"^@event.pt.type,
"ENREGRRCR" - 001^@event.oth.vendor.tx.type,"Y"^@event.undoable,"RCR"^@event.pt.type,
"ENTER REG RCR" ,"eereg"^@event.stats.program,
"EDREGRRCR" - 004^@event.oth.vendor.tx.type,"RCR"^@event.pt.type,
"EDIT REG RCR" ,"eereg"^@event.stats.program,
"ENDISRRCR" - "DISCHARGE RCR" ,003^@event.oth.vendor.tx.type,
"Y"^@event.undoable,"RCR"^@event.pt.type,"endisrcr"^@event.stats.program,
"UNDODISRRCR" - "UNDO DISCHARGE RCR" ,
012^@event.oth.vendor.tx.type,"RCR"^@event.pt.type,"undodisrcr"^@event.stats.program,
"EDDISRRCR" - 004^@event.oth.vendor.tx.type,"RCR"^@event.pt.type,
"EDIT DIS RCR" ,"eereg"^@event.stats.program,
"FXREGRRCR" - 006^@event.oth.vendor.tx.type,"RCR"^@event.pt.type,
"UNDO RCR REGISTRATION" ,"fxreg"^@event.stats.program,
"REVRRCR" - 014^@event.oth.vendor.tx.type,"RCR"^@event.pt.type,
"ADD RCR REVISIT" ,"revrcr"^@event.stats.program,
"DELREVRRCR" - 015^@event.oth.vendor.tx.type,"RCR"^@event.pt.type,
"DELETE RCR REVISIT" ,"revrcr"^@event.stats.program,
"CANPRERCR" - "CANCEL PRE RCR" ,006^@event.oth.vendor.tx.type,"RCR"^@event.pt.type,
"CANSCHRRCR" - "CANCEL SCH RCR" ,006^@event.oth.vendor.tx.type,"RCR"^@event.pt.type,

RCR-Outpatient room/bed events

"ENRCRBED" - 002^@event.oth.vendor.tx.type,"RCR"^@event.pt.type,"Y"^@event.undoable,
"PLACE RCR IN ROOM/BED" ,"enrcrbed"^@event.stats.program,
"EDRCRBED" - 002^@event.oth.vendor.tx.type,"RCR"^@event.pt.type,"Y"^@event.undoable,
"EDIT RCR ROOM/BED" ,"edrcrbed"^@event.stats.program,
"RMRCRBED" - 002^@event.oth.vendor.tx.type,"RCR"^@event.pt.type,"Y"^@event.undoable,
"REMOVE RCR FROM ROOM/BED" ,"rmrcrbed"^@event.stats.program,
"FXRCRRRACCOM" - ""^@event.oth.vendor.tx.type,"RCR"^@event.pt.type,"Y"^@event.undoable,
"EDIT RCR ROOM RATE ACCOM" ,"fxrcrrraccom"^@event.stats.program,
"FXRCRSVC" - ""^@event.oth.vendor.tx.type,"RCR"^@event.pt.type,"Y"^@event.undoable,
"EDIT RCR SERVICE" ,"fxcrsvc"^@event.stats.program,
"FXRCRACCOM" - ""^@event.oth.vendor.tx.type,"RCR"^@event.pt.type,"Y"^@event.undoable,
"EDIT RCR REQUESTED ACCOM" ,"fxrcraccom"^@event.stats.program,
"UNDOENRCRBED" - 002^@event.oth.vendor.tx.type,"RCR"^@event.pt.type,
"UNDO PLACE RCR IN ROOM/BED" ,"undoenrcrbed"^@event.stats.program,
"UNDOEDRCRBED" - 002^@event.oth.vendor.tx.type,"RCR"^@event.pt.type,
"UNDO EDIT RCR ROOM/BED" ,"undoedopbed"^@event.stats.program,
"UNDOFXRCRSVC" - ""^@event.oth.vendor.tx.type,"RCR"^@event.pt.type,

"UNDO EDIT RCR SERVICE" , "undofxrcrsvc"^@event.stats.program,
"UNDOFXRCRACCOM" - ""^@event.oth.vendor.tx.type,"RCR"^@event.pt.type,
"UNDO EDIT RCR REQUESTED ACCOM" , "undofxrcraccom"^@event.stats.program,
"UNDOFXRCRRRACCOM" - ""^@event.oth.vendor.tx.type,"RCR"^@event.pt.type,
"UNDO EDIT RCR ROOM RATE ACCOM" , "undofxrcrrraccom"^@event.stats.program,
"UNDORMRCRBED" - ""^@event.oth.vendor.tx.type,"RCR"^@event.pt.type,
"UNDO REMOVE RCR FROM ROOM/BED" , "undormrcrbed"^@event.stats.program,

REF outpatients

"ENSCHREF" - "ENTER SCH REF" ,001^@event.oth.vendor.tx.type,"REF"^@event.pt.type,
"EDSCHREF" - "EDIT SCH REF" ,004^@event.oth.vendor.tx.type,"REF"^@event.pt.type,
"ENPREREF" - "ENTER PRE REF" ,001^@event.oth.vendor.tx.type,"REF"^@event.pt.type,
"EDPREREF" - "EDIT PRE REF" ,004^@event.oth.vendor.tx.type,"REF"^@event.pt.type,
"ENREGREF" - 001^@event.oth.vendor.tx.type,"Y"^@event.undoable,"REF"^@event.pt.type,
"ENTER REG REF" , "eereg"^@event.stats.program,
"EDREGREF" - 004^@event.oth.vendor.tx.type,"REF"^@event.pt.type,
"EDIT REG REF" , "eereg"^@event.stats.program,
"FXREGREF" - 006^@event.oth.vendor.tx.type,"REF"^@event.pt.type,
"UNDO REF REGISTRATION" , "fxreg"^@event.stats.program,
"CANPREREF" - "CANCEL PRE REF" ,006^@event.oth.vendor.tx.type,"REF"^@event.pt.type,
"CANSCHREF" - "CANCEL SCH REF" ,006^@event.oth.vendor.tx.type,"REF"^@event.pt.type,

REF - Outpatient room/bed events

"ENREFBED" - 002^@event.oth.vendor.tx.type,"REF"^@event.pt.type,"Y"^@event.undoable,
"PLACE REF IN ROOM/BED" , "enopbed"^@event.stats.program,
"EDREFBED" - 002^@event.oth.vendor.tx.type,"REF"^@event.pt.type,"Y"^@event.undoable,
"EDIT REF ROOM/BED" , "edopbed"^@event.stats.program,
"RMREFBED" - 002^@event.oth.vendor.tx.type,"REF"^@event.pt.type,"Y"^@event.undoable,
"REMOVE REF FROM ROOM/BED" , "rmopbed"^@event.stats.program,
"ENDEPREF" - 003^@event.oth.vendor.tx.type,"REF"^@event.pt.type,"Y"^@event.undoable,
"ENTER REF DEPARTURE" , "endepout"^@event.stats.program,
"EDDEPREF" - 004^@event.oth.vendor.tx.type,"REF"^@event.pt.type,
"EDIT DEPARTED REF" , "eereg"^@event.stats.program,
"FXDEPREFDFT" - ""^@event.oth.vendor.tx.type,"REF"^@event.pt.type,
"EDIT REF DEPART DATE" , "fxdepdt"^@event.stats.program,
"FXREFRRACCOM" - ""^@event.oth.vendor.tx.type,"REF"^@event.pt.type,"Y"^@event.undoable,
"EDIT REF ROOM RATE ACCOM" , "fxoprpraccom"^@event.stats.program,
"FXREFSVC" - ""^@event.oth.vendor.tx.type,"REF"^@event.pt.type,"Y"^@event.undoable,
"EDIT REF SERVICE" , "fxopsvc"^@event.stats.program,
"FXREFACCOM" - ""^@event.oth.vendor.tx.type,"REF"^@event.pt.type,"Y"^@event.undoable,
"EDIT REF REQUESTED ACCOM" , "fxopaccom"^@event.stats.program,
"UNDOENREFBED" - 002^@event.oth.vendor.tx.type,"REF"^@event.pt.type,

"UNDO PLACE REF IN ROOM/BED" ,"undoenopbed"^@event.stats.program,
"UNDOEDREFBED" - 002^@event.oth.vendor.tx.type,"REF"^@event.pt.type,
"UNDO EDIT REF ROOM/BED" ,"undoedopbed"^@event.stats.program,
"UNDOFXREFSVC" - ""^@event.oth.vendor.tx.type,"REF"^@event.pt.type,
"UNDO EDIT REF SERVICE" ,"undofxopsvc"^@event.stats.program,
"UNDOFXREFACCOM" - ""^@event.oth.vendor.tx.type,"REF"^@event.pt.type,
"UNDO EDIT REF REQUESTED ACCOM" ,"undofxopaccom"^@event.stats.program,
"UNDOFXREFRRACCOM" - ""^@event.oth.vendor.tx.type,"REF"^@event.pt.type,
"UNDO EDIT REF ROOM RATE ACCOM" ,"undofxopraccom"^@event.stats.program,
"UNDORMREFBED" - ""^@event.oth.vendor.tx.type,"REF"^@event.pt.type,
"UNDO REMOVE REF FROM ROOM/BED" ,"undormopbed"^@event.stats.program,
"UNDODEPREF" - 012^@event.oth.vendor.tx.type,"REF"^@event.pt.type,
"UNDO REF DEPARTURE" ,"undodepout"^@event.stats.program,

SDC outpatients

"ENSCHSDC" - "ENTER SCH SDC" ,001^@event.oth.vendor.tx.type,"SDC"^@event.pt.type,
"EDSCHSDC" - "EDIT SCH SDC" ,004^@event.oth.vendor.tx.type,"SDC"^@event.pt.type,
"ENPRESDC" - "ENTER PRE SDC" ,001^@event.oth.vendor.tx.type,"SDC"^@event.pt.type,
"EDPRESDC" - "EDIT PRE SDC" ,004^@event.oth.vendor.tx.type,"SDC"^@event.pt.type,
"ENREGSDC" - 001^@event.oth.vendor.tx.type,"Y"^@event.undoable,"SDC"^@event.pt.type,
"ENTER REG SDC" ,"eereg"^@event.stats.program,
"EDREGSDC" - 004^@event.oth.vendor.tx.type,"SDC"^@event.pt.type,
"EDIT REG SDC" ,"eereg"^@event.stats.program,
"FXREGSDC" - 006^@event.oth.vendor.tx.type,"SDC"^@event.pt.type,
"UNDO SDC REGISTRATION" ,"fxreg"^@event.stats.program,
"CANPRESDC" - "CANCEL PRE SDC" ,006^@event.oth.vendor.tx.type,"SDC"^@event.pt.type,
"CANSCHSDC" - "CANCEL SCH SDC" ,006^@event.oth.vendor.tx.type,"SDC"^@event.pt.type,

SDC-Outpatient room/bed events

"ENSDCBED" - 002^@event.oth.vendor.tx.type,"SDC"^@event.pt.type,"Y"^@event.undoable,
"PLACE SDC IN ROOM/BED" ,"enopbed"^@event.stats.program,
"EDSDCBED" - 002^@event.oth.vendor.tx.type,"SDC"^@event.pt.type,"Y"^@event.undoable,
"EDIT SDC ROOM/BED" ,"edopbed"^@event.stats.program,
"RMSDCBED" - 002^@event.oth.vendor.tx.type,"SDC"^@event.pt.type,"Y"^@event.undoable,
"REMOVE SDC FROM ROOM/BED" ,"rmopbed"^@event.stats.program,
"ENDEPSDC" - 003^@event.oth.vendor.tx.type,"SDC"^@event.pt.type,"Y"^@event.undoable,
"ENTER SDC DEPARTURE" ,"endepout"^@event.stats.program,
"EDDEPSDC" - 004^@event.oth.vendor.tx.type,"SDC"^@event.pt.type,
"EDIT DEPARTED SDC" ,"eereg"^@event.stats.program,
"FXDEPSDCDT" - ""^@event.oth.vendor.tx.type,"SDC"^@event.pt.type,
"EDIT SDC DEPART DATE" ,"fxdepdt"^@event.stats.program,
"FXSDCRRACCOM" - ""^@event.oth.vendor.tx.type,"SDC"^@event.pt.type,"Y"^@event.undoable,

"EDIT SDC ROOM RATE ACCOM" ,"fxoprracom"^@event.stats.program,
"FXSDCSVC" - ""^@event.oth.vendor.tx.type,"SDC"^@event.pt.type,"Y"^@event.undoable,
"EDIT SDC SERVICE" ,"fxopsvc"^@event.stats.program,
"FXSDCACCOM" - ""^@event.oth.vendor.tx.type,"SDC"^@event.pt.type,"Y"^@event.undoable,
"EDIT SDC REQUESTED ACCOM" ,"fxopacom"^@event.stats.program,
"UNDOENSDCBED" - 002^@event.oth.vendor.tx.type,"SDC"^@event.pt.type,
"UNDO PLACE SDC IN ROOM/BED" ,"undoenopbed"^@event.stats.program,
"UNDOEDSDCBED" - 002^@event.oth.vendor.tx.type,"SDC"^@event.pt.type,
"UNDO EDIT SDC ROOM/BED" ,"undodopbed"^@event.stats.program,
"UNDOFXSDCSVC" - ""^@event.oth.vendor.tx.type,"SDC"^@event.pt.type,
"UNDO EDIT SDC SERVICE" ,"undofxopsvc"^@event.stats.program,
"UNDOFXSDCACCOM" - ""^@event.oth.vendor.tx.type,"SDC"^@event.pt.type,
"UNDO EDIT SDC REQUESTED ACCOM" ,"undofxopacom"^@event.stats.program,
"UNDOFXSDCRRACCOM" - ""^@event.oth.vendor.tx.type,"SDC"^@event.pt.type,
"UNDO EDIT SDC ROOM RATE ACCOM" ,"undofxoprracom"^@event.stats.program,
"UNDORMSDCBED" - ""^@event.oth.vendor.tx.type,"SDC"^@event.pt.type,
"UNDO REMOVE SDC FROM ROOM/BED" ,"undormopbed"^@event.stats.program,
"UNDODEPSDC" - 012^@event.oth.vendor.tx.type,"SDC"^@event.pt.type,
"UNDO SDC DEPARTURE" ,"undodepout"^@event.stats.program,

POV outpatients

"ENSCHPOV" - "ENTER SCH POV" ,001^@event.oth.vendor.tx.type,"POV"^@event.pt.type,
"EDSCHPOV" - "EDIT SCH POV" ,004^@event.oth.vendor.tx.type,"POV"^@event.pt.type,
"ENPREPOV" - "ENTER PRE POV" ,001^@event.oth.vendor.tx.type,"POV"^@event.pt.type,
"EDPREPOV" - "EDIT PRE POV" ,004^@event.oth.vendor.tx.type,"POV"^@event.pt.type,
"ENREGPOV" - 001^@event.oth.vendor.tx.type,"Y"^@event.undoable,"POV"^@event.pt.type,
"ENTER REG POV" ,"eereg"^@event.stats.program,
"EDREGPOV" - 004^@event.oth.vendor.tx.type,"POV"^@event.pt.type,
"EDIT REG POV" ,"eereg"^@event.stats.program,
"FXREGPOV" - 006^@event.oth.vendor.tx.type,"POV"^@event.pt.type,
"UNDO POV REGISTRATION" ,"fxreg"^@event.stats.program,
"CANPREPOV" - "CANCEL PRE POV" ,006^@event.oth.vendor.tx.type,"POV"^@event.pt.type,
"CANSCHPOV" - "CANCEL SCH POV" ,006^@event.oth.vendor.tx.type,"POV"^@event.pt.type,

POV-Outpatient room/bed events

"ENPOVBED" - 002^@event.oth.vendor.tx.type,"POV"^@event.pt.type,"Y"^@event.undoable,
"PLACE POV IN ROOM/BED" ,"enopbed"^@event.stats.program,
"EDPOVBED" - 002^@event.oth.vendor.tx.type,"POV"^@event.pt.type,"Y"^@event.undoable,
"EDIT POV ROOM/BED" ,"edopbed"^@event.stats.program,
"RMPOVBED" - 002^@event.oth.vendor.tx.type,"POV"^@event.pt.type,"Y"^@event.undoable,
"REMOVE POV FROM ROOM/BED" ,"rmopbed"^@event.stats.program,
"ENDEPPOV" - 003^@event.oth.vendor.tx.type,"POV"^@event.pt.type,"Y"^@event.undoable,

"ENTER POV DEPARTURE" ,"endepout"^@event.stats.program,
"EDDEPPOV" - 004^@event.oth.vendor.tx.type,"POV"^@event.pt.type,
"EDIT DEPARTED POV" ,"eereg"^@event.stats.program,
"FXDEPPOVDT" - ""^@event.oth.vendor.tx.type,"POV"^@event.pt.type,
"EDIT POV DEPART DATE" ,"fxdepdt"^@event.stats.program,
"FXPOVRRACCOM" - ""^@event.oth.vendor.tx.type,"POV"^@event.pt.type,"Y"^@event.undoable,
"EDIT POV ROOM RATE ACCOM" ,"fxopraccom"^@event.stats.program,
"FXPOVSV" - ""^@event.oth.vendor.tx.type,"POV"^@event.pt.type,"Y"^@event.undoable,
"EDIT POV SERVICE" ,"fxopsv" "^@event.stats.program,
"FXPOVACCOM" - ""^@event.oth.vendor.tx.type,"POV"^@event.pt.type,"Y"^@event.undoable,
"EDIT POV REQUESTED ACCOM" ,"fxopaccom"^@event.stats.program,
"UNDOENPOVBED" - 002^@event.oth.vendor.tx.type,"POV"^@event.pt.type,
"UNDO PLACE POV IN ROOM/BED" ,"undoenopbed"^@event.stats.program,
"UNDOEDPOVBED" - 002^@event.oth.vendor.tx.type,"POV"^@event.pt.type,
"UNDO EDIT POV ROOM/BED" ,"undodopbed"^@event.stats.program,
"UNDOFXPOVSV" - ""^@event.oth.vendor.tx.type,"POV"^@event.pt.type,
"UNDO EDIT POV SERVICE" ,"undofxopsv" "^@event.stats.program,
"UNDOFXPOVACCOM" - ""^@event.oth.vendor.tx.type,"POV"^@event.pt.type,
"UNDO EDIT POV REQUESTED ACCOM" ,"undofxopaccom"^@event.stats.program,
"UNDOFXPOVRRACCOM" - ""^@event.oth.vendor.tx.type,"POV"^@event.pt.type,
"UNDO EDIT POV ROOM RATE ACCOM" ,"undofxopraccom"^@event.stats.program,
"UNDORMPOVBED" - ""^@event.oth.vendor.tx.type,"POV"^@event.pt.type,
"UNDO REMOVE POV FROM ROOM/BED" ,"undormopbed"^@event.stats.program,
"UNDODEPPOV" - 012^@event.oth.vendor.tx.type,"POV"^@event.pt.type,
"UNDO POV DEPARTURE" ,"undodepout"^@event.stats.program,

PPR outpatients

"ENSCHPPR" - "ENTER SCH PPR" ,""^@event.oth.vendor.tx.type,"PPR"^@event.pt.type,
"EDSCHPPR" - "EDIT SCH PPR" ,""^@event.oth.vendor.tx.type,"PPR"^@event.pt.type,
"ENPREPPR" - "ENTER PRE PPR" ,""^@event.oth.vendor.tx.type,"PPR"^@event.pt.type,
"EDPREPPR" - "EDIT PRE PPR" ,""^@event.oth.vendor.tx.type,"PPR"^@event.pt.type,
"ENREGPPR" - ""^@event.oth.vendor.tx.type,"Y"^@event.undoable,"PPR"^@event.pt.type,"ENTER
REG PPR" ,
"EDREGPPR" - ""^@event.oth.vendor.tx.type,"PPR"^@event.pt.type,"EDIT REG PPR" ,
"FXREGPPR" - ""^@event.oth.vendor.tx.type,"PPR"^@event.pt.type,"UNDO PPR REGISTRATION" ,
"CANPREPPR" - "CANCEL PRE PPR" ,""^@event.oth.vendor.tx.type,"PPR"^@event.pt.type,
"CANSCHPPR" - "CANCEL SCH PPR" ,
"ENDISPPR" - ""^@event.oth.vendor.tx.type,"PPR"^@event.pt.type,"Y"^@event.undoable,"ENTER
PPR DISCHARGE" ,
"EDDISPPR" - ""^@event.oth.vendor.tx.type,"PPR"^@event.pt.type,"EDIT DISCHARGED PPR" ,
"FXDISPPRDT" - ""^@event.oth.vendor.tx.type,"PPR"^@event.pt.type,"EDIT PPR DISCHARGE
DATE" ,

"UNDODISPPR" - ""^@event.oth.vendor.tx.type,"PPR"^@event.pt.type,"UNDO PPR DISCHARGE" ,
"ENCPPR" - 001^@event.oth.vendor.tx.type,"PPR"^@event.pt.type,"ENTER PPR ENCOUNTER" ,
"AUTHPPR" - ""^@event.oth.vendor.tx.type,"PPR"^@event.pt.type,"ENTER PPR AUTHORIZATION" ,
END;

About D288 IT Solutions



Our mission at D288 IT Solutions is to diligently work alongside healthcare organizations for the advancement of integration.

Check us out at <http://www.d288itsolutions.com>

Connect with us at

Facebook: www.facebook.com/d288itsolutions

Twitter: www.twitter.com/d288itsolutons

LinkedIn: www.linkedin.com/companies/d288-it-solutions

Email: info@d288itsolutions.com

